

# Galaxy for virologist training

## Exercise 4: Nanopore mapping and Assembly 101

Title	Galaxy
<b>Training dataset:</b>	Nanopore MinION Sequencing of a Monkey Pox Virus (MPXV) from Spain 2022 outbreak. Data is publicly available at SRA with <a href="#">ID ERR10297654</a> . <a href="#">Paper</a>
<b>Questions:</b>	<ul style="list-style-type: none"><li>• How Nanopore reads are differently assembled from Illumina?</li></ul>
<b>Objectives:</b>	<ul style="list-style-type: none"><li>• Understand the concept of assembly</li><li>• Learn how to interpret assembly quality control metrics</li></ul>
<b>Estimated time:</b>	40 min

## 1. Description


Nanopore technology is a third generation sequencing technique which allows to get longer sequences, but with reduced sequence quality. Different technologies have different formats, qualities, and specific known biases which make the analysis different among them. In this tutorial, we are going to see an example of how to assemble long reads from a Nanopore sequencing run.

## 2. Upload data to galaxy

### Training dataset

- [SRA ID: ERR10297654]([https://trace.ncbi.nlm.nih.gov/Traces/?view=run\\_browser&acc=ERR10297654&display=metadata](https://trace.ncbi.nlm.nih.gov/Traces/?view=run_browser&acc=ERR10297654&display=metadata))

### Create new history

- Click the  icon at the top of the history panel and create a new history with the name `nanopore assembly 101 tutorial` as explained [here](#)

### Upload data

1. Look for `SRA` in the tool search bar and select `Faster Download and Extract Reads in FASTQ format from NCBI SRA`
2. Accession = `ERR10297654`
3. Execute

## Table of Contents

- [1. Description](#)
- [2. Upload data to galaxy](#)
  - [Training dataset](#)
  - [Create new history](#)
  - [Upload data](#)
  - [Load reference file from NCBI](#)
  - [Unhide data](#)
  - [Mapping with Minimap2](#)
  - [Mapping stats with samtools](#)
  - [Assemble reads with Flye](#)
  - [Assembly quality control with Quast](#)

**Galaxy Europe** Workflow Visualize Shared Data Help User

**Tools**

SRA 1

Upload Data

Show Sections

SRA server

Download and Generate Pileup

Format from NCBI SRA

Download and Extract Reads in BAM format from NCBI SRA

**Faster Download and Extract Reads in FASTQ format from NCBI SRA** 2

pysradb search sequence metadata from SRA/ENA

EBI SRA ENA SRA

LowMemPeakPickerHiResRandomAccess Finds mass spectrometric peaks in profile mass spectra.

Make.sra creates the necessary files for a NCBI submission

**Faster Download and Extract Reads in FASTQ format from NCBI SRA (Galaxy Version 2.11.0+galaxy1)**

select input type

SRR accession

**Accession**

ERR10297654

Must start with SRR, DRR or ERR, e.g. SRR925743, ERR343809

Advanced Options

Email notification

No

Send an email notification when the job completes.

Execute 4

**What it does?**

This tool extracts data (in fastq format) from the Short Read Archive (SRA) at the National Center for Biotechnology Information (NCBI). It is based on the fasterq-dump utility of the SRA Toolkit.

**How to use it?**

## Load reference file from NCBI

1. Search NCBI using the search toolbox and select NCBI Accession Download Download sequences from GenBank/RefSeq by accession through the NCBI ENTREZ API
2. Select source for IDs > Direct entry
3. ID List = NC\_063383.1
4. Execute

**Galaxy Europe** Flujo de Trabajo Visualizar Datos Compartidos Ayuda Usuario

Herramientas

ncbi 1

Cargar Datos

Show Sections

Krona pie chart from taxonomic profile

Kraken taxonomic report view report of classification for multiple samples

Unipept retrieve taxonomy for peptides

MaxQuant (using mqpar.xml)

Krona pie chart from taxonomic profile

NCBI ESearch search NCBI Databases by text query

NCBI EPost post UIDs to NCBI History Server

NCBI EFetch fetch records from NCBI

**NCBI Accession Download Download sequences from GenBank/RefSeq by accession through the NCBI ENTREZ API** 2

NCBI EInfo fetch NCBI database metadata

NCBI ECitMatch search NCBI for citations in PubMed

ETE tax DB parameter generate the

**NCBI Accession Download Download sequences from GenBank/RefSeq by accession through the NCBI ENTREZ API (Galaxy Version 0.2.7+galaxy0)**

Select source for IDs

Direct Entry 3

ID List

NC\_063383.1 4

Newline/Comma separated list of IDs

Molecule Type

Nucleotide

File Format

FASTA

How to handle download failures

Abort with error on first failure

Add accession to failed list and continue

Email notification

No

Send an email notification when the job completes.

Execute 5

History

buscar conjuntos de datos

nanopore assembly 101 tutorial

361 MB

4 : fasterq-dump log

3 : Other data (fasterq-dump)

a list with 0 datasets

2 : Single-end data (fasterq-dump)

a list with 1 fastqsanger.gz dataset

1 : Pair-end data (fasterq-dump)

a list with 0 pairs

## Unhide data

Using SRA and NCBI API downloads data as hidden so we are going to unhide this data as follows: 1. Click on the strikethrough eye (Show hidden) 2. Select the strikethrough for ERR10297654 and NC\_063383.1 datas. 3. Then select the location icon (show active)


## Mapping with Minimap2

1. Search `minimap2` using the search toolbox and select `Map with minimap2 A fast pairwise aligner for genomic and spliced nucleotide sequences`
2. Will you select a reference genome from your history or use a built-in index?: Use a genome from history and built-in index
  - Select NC\_063383.1
3. Select fastq dataset: ERR10297654
4. Select a profile of preset options > Oxford Nanopore Read to reference mapping (map-ont)
5. Click execute and wait.

## Mapping stats with samtools

1. Search `flagstat` using the search toolbox and select `Samtools flagstat tabulate descriptive stats for BAM dataset`
2. BAM File to report statistics of > Select Minimap2 bam output
3. Click execute and wait.

4. Click in the   and see the bam stats.



**Herramientas**

flagstats 1

Cargar Datos

Show Sections

flagstat provides simple stats on BAM files

Samtools flagstat tabulate descriptive stats for BAM dataset 2

**FLUJOS DE TRABAJO**

Todos los flujos de trabajo

**Samtools flagstat tabulate descriptive stats for BAM dataset (Galaxy Version 2.0.4)**

BAM File to report statistics of

9: Map with minimap2 on data 5 and data 8 (mapped reads in BAM format) 3

Output format

txt

(--output-fmt)

Email notification

No


Send an email notification when the job completes.

Execute 4

- Which is the mapping rate?
- How many reads do we have in our dataset?

## Assemble reads with Flye

1. Search Flye assembler using the search toolbox and select Flye de novo assembler for single molecule sequencing reads
2. Input reads: ERR10297654
3. Click execute and wait.



**Herramientas**

Flye 1

Cargar Datos

Show Sections

Flye de novo assembler for single molecule sequencing reads 2

**FLUJOS DE TRABAJO**

Todos los flujos de trabajo

**Flye de novo assembler for single molecule sequencing reads (Galaxy Version 2.9.1+galaxy0)**

Input reads

8: NC\_063383.1  
5: ERR10297654 3

Mode

Nanopore raw (--nano-raw)

Number of polishing iterations

1

Polishing is performed as the final assembly stage. By default, Flye runs one polishing iteration. Additional iterations might correct a small number of extra errors (due to improvements on how reads may align to the corrected assembly). If the parameter is set to 0, the polishing is not performed (--iterations)

Minimum overlap between reads

4 Execute

## Assembly quality control with Quast

1. Search Quast in the search tool box.
2. Contigs/scaffolds file > Flye results
3. Use a reference genome: Yes. Select the NC\_063383.1 fasta file previously loaded.

**Galaxy Europe** Flujo de Trabajo Visualizar Datos Compartidos Ayuda Usuario

**Herramientas**

quast 1

Cargar Datos

Show Sections

**Quast Genome assembly Quality** 2

rnaQUAST A Quality Assessment Tool for De Novo Transcriptome Assemblies

FLUJOS DE TRABAJO

Todos los flujos de trabajo

**Quast Genome assembly Quality (Galaxy Version 5.2.0+galaxy0)**

Use customized names for the input files?

No, use dataset names

They will be used in reports, plots and logs

**Contigs/scaffolds file**

11: Flye on data 5: consensus  
8: NC\_063383.1

**Reads options**

Disabled

Currently, the supported read types are Illumina unpaired, paired-end and mate-pair reads, PacBio SMRT, and Oxford Nanopore long reads.

**Type of assembly**

Genome

**Use a reference genome?**

Yes

Many metrics can't be evaluated without a reference. If this is omitted, QUAST will only report the metrics that can be evaluated without a reference.

**Reference genome**

11: Flye on data 5: consensus  
8: NC\_063383.1

**6 Execute**

1. Click the :eye: icon Quast HTML report.

- How much of or reference genome have we reconstructed?
- How many contigs do we have greater than 1000 pb?
- How long is the largest contig in the assembly?
- Which is the N50?

2. Open the Icarus viewer in the quast report.

**QUAST**  
Quality Assessment Tool for Genome Assemblies by CAB

17 November 2021, Wednesday, 18:07:45

[View in Icarus contig browser](#)

All statistics are based on contigs of size  $\geq 500$  bp, unless otherwise noted (e.g., "# cor bp)" and "Total length ( $\geq 0$  bp)" include all contigs).

Aligned to "dataset\_76f460df\_9dce\_4919\_9108\_be70c4d29af9" | 11 029 bp | 1 fragment G+C

**Genome statistics** SPAdes\_on\_data\_3\_and\_data\_2\_...

Genome fraction (%) 83.070

► Which contig align against our reference genome?

This training history is available at:

<https://usegalaxy.eu/u/s.varona/h/nanopore-assembly-101-tutorial>

Note: Nanopore data is known to have more error than short sequencing reads. This is why assembly post-processing is strongly recommended, usually using combined sequencing approximation with both Nanopore and Illumina reads.